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
Pro Cellular Body Scale      BS-2001-G1

CE

- Thank you very much for selecting our Pro Cellular Body Scale.
- Please read the entire user manual carefully to ensure the safe usage of this product, and keep the manual for future reference in case you have questions.

# Table of Contents

- Symbols** ..... 2
- Safety Information**
- Where to Use and Store Your Scale ..... 2
- Proper Use of Your Scale ..... 2
- Caution ..... 3
- Overview**
- Device Components ..... 4
- List of Parts ..... 4
- Scale Display Codes ..... 5
- Getting Started**
- Installing the Batteries and Pairing the Scale ..... 6
- Select Unit of Measure ..... 7
- Start Measuring** ..... 8
- Troubleshooting**
- Error Prompt ..... 9
- When Measuring ..... 10
- Specifications** ..... 11
- Maintenance** ..... 12
- Warranty** ..... 12
- Contact Information** ..... 13

♥ Symbols

	Manufacturer		Date of manufacture
<b>SN</b>	Serial number		
	Direct current		Do not dispose of this product as unsorted municipal waste, take it for recycling.
	Type BF applied part (top surface)		
	Caution: these notes must be observed to prevent any damage to the device.		Recycle
<b>CE</b>	Symbol for "COMPLIES WITH EU REQUIREMENTS"		

♥ Where to Use and Store Your Scale

To ensure your safety as well as the service life of your scale, avoid using the scale under the following conditions:

- Slippery floor such as tile floor
- Jumping onto the platform immediately after bath or with wet feet
- Near a cell phone or microwave oven

Avoid storing your scale in the following locations:

- Where there is water
- Where the device may be exposed to extreme temperatures, humidity, moisture, direct sunlight, or dust
- Where there is risk of the scale being shocked or dropped
- Near where chemicals or corrosive gases are stored
- Within the reach of infants or children

♥ Proper Use of Your Scale


To ensure the accuracy of measurements, please follow the instructions below before you weigh yourself.

- Place the scale on a hard, flat surface. Soft surfaces such as carpet will affect the scale's accuracy.
- Weigh yourself at least two hours after eating, preferably after getting up in the morning.
- Avoid weighing yourself immediately after strenuous exercise, taking a sauna or bath, drinking, or eating.
- Always weigh yourself at the same time of the day, and on the same scale located on the same hard, flat surface.

 CAUTION

- This device is designed for use on people 10 years and older and is not intended for infants.
- This device is designed for use in the home healthcare environment to measure a person's weight.
- This device is suitable for people who can read and understand this user guide. The device should not be used to diagnose or treat any medical condition. Any information provided by this device is not meant to treat, cure, or prevent any disease or illness from happening. If in doubt, contact your doctor.
- The scale is equipped with a data transmission function. It may emit electromagnetic energy so as to perform its intended function. Nearby portable and mobile RF communications equipment can affect the performance of the scale. And the scale can interfere with electrical equipment in its vicinity.
- Please use the device according to the user manual. Any misuse can cause electric shock, burns, fire, and other unexpected hazards.
- Please use and store the device per the environmental conditions provided in the user manual.
- Don't expose the device to extreme temperatures, direct sunlight, moist, or corrosive environment.
- Don't step on the scale when your body or feet are wet, especially after bathing or showering, to prevent slipping.
- Please keep the device out of reach of infants, children, and pets, since inhalation or swallowing of small parts is dangerous or even fatal.
- Don't use the device beyond the expected lifetime because degraded sensors can degrade performance or cause other problems.
- The patient is an intended operator. The patient can measure, transmit data, and change batteries under normal circumstances and maintain the device and its accessories according to the user manual.
- The scale requires 30 minutes to warm from the minimum and cool from maximum storage temperature between uses until it is ready for measurement.
- After drop or shock that may cause changes in the performance, please contact Customer Support.
- Don't open or repair the device by yourself. The use of heat- and cold-producing devices, such as electric heating blankets, heating pads, or ice packs, may impair the performance of device and increase the risk of injury to the user.
- If you are allergic to plastic or rubber, don't use this device.
- Do not use the device if it is damaged in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.
- WARNING: Don't use this device in an oxygen-rich environment.
- WARNING: Don't use near active High Frequency surgical equipment and the RF shielded room of an Medical Equipment system for magnetic resonance imaging, where the intensity of Electro Magnetic disturbances is high.
- WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of this equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- WARNING: No modifications of this equipment are allowed. This may result in increased EMISSIONS or decreased IMMUNITY of the scale.


♥ Device Components

♥ List of Parts

1. Pro Cellular Body Scale
2. Four AA-size Batteries (1.5V each)
3. User Manual

♥ Scale Display Codes


<b>kg</b>	Kilogram	✕	Signal symbol
<b>st lb</b>	Stone		Data transmitting
<b>lb</b>	Pound		Battery symbol


## 1. Installing the Batteries and Pairing the Scale

- Open the battery door on the back of the scale.
- Insert the batteries (4 x 1.5V AA) into the battery compartment according to the polarity indications marked inside the compartment.

- The scale will automatically pair-up with the system when used for the first time. The signal symbol will be displayed at the same time. It indicates the quality of the current network signal.

If the device successfully pairs up, the symbol “[ ]” will be shown on the LCD.

If the device fails to pair up, the symbol “E1” will be shown on the LCD.

### ⚠ CAUTION

- When the symbol “” and “Lo” appear, the device will power off. Then replace with a new set of batteries. Please replace all four batteries at the same time. Do NOT mix the old batteries with the new one.
- Worn batteries are hazardous waste. Do NOT dispose of them together with the household garbage.
- Please refer to the local ordinances and recycling instructions regarding disposal of the worn batteries and scrapped device.
- Do not mix alkaline, standard (carbon-zinc), or rechargeable (Ni-Cd, Ni-MH, etc) batteries.
- If you do not intend to use this unit for a prolonged period of time, it is advisable to remove the batteries before storing.

-- CAUTION! RISK OF EXPLOSION IF BATTERY IS REPLACED BY AN INCORRECT TYPE.

-- DISPOSE OF USED BATTERIES ACCORDING TO THE INSTRUCTIONS.


-- Disposal of a battery into fire or a hot oven, or mechanically crushing or cutting of a battery, that can result in an explosion.

-- Leaving a battery in an extremely high temperature surrounding environment that can result in an explosion or the leakage of flammable liquid or gas;

-- A battery subjected to extremely low air pressure that may result in an explosion or the leakage of flammable liquid or gas.


## 2. Select Unit of Measure

- With batteries correctly installed, you will see the black “UNIT” button into the battery door. Press it to select your preferred unit of measurement. The default measurement unit is “lb” You may press the “UNIT” button to choose between pounds, kilograms, and stones.

- Close battery door

1. Step on the platform of the device and remain still.


2. Do not move or talk during the measurement.

3. When the measurement concludes, the display briefly blinks and results are shown on the display.

4. Step off the scale platform.

5. The device automatically transmits measurement results.  
When data is successfully transmitted, "OK" will be displayed and icon "↻" stops blinking.

## ♥ Error Prompt

Code	Description	Resolution
	Overload	This device is only accurate up to its rated capacity.
	Battery low	1) Replace batteries with four Duracell AAs. 2) Check that batteries are all brand new and not expired.
	Cellular connection error or data transmission error	1) Move the device to another area, preferably closer to a window. 2) If the new area has strong cellular signal and the issues persists, contact customer support.
	Signal not strong enough	1) Replace batteries with four Duracell AAs. 2) Move the device to a strong signal place, preferably closer to a window.
	SIM error	Contact customer support.
	IC is damaged	Contact customer support.

♥ When Measuring

Problem	Possible reason	Resolution
No display	Batteries not yet installed or batteries are not installed correctly.	Install the batteries. (Please refer to Installing the Batteries and Pairing the Scale.)
	Battery is running low.	Replace batteries with four Duracell AAs.
The device powers off automatically.	Low battery.	Replace batteries with four Duracell AAs.

♥ Specifications

<b>Product Name</b>	Pro Cellular Body Scale	
<b>Dimension</b>	Scale: Approx. 360x360x37.2mm	
<b>Net Weight</b>	Approximately 6.7lb ( Excluding the dry cells)	
<b>Display</b>	Digital LCD	
<b>Measurement Unit</b>	Pound / Kilogram / Stone	
<b>Measurement Range</b>	11lb to 551.2lb / 5kg to 250kg / 0st:11lb to 39st: 5lb	
<b>Division</b>	0.1kg / 0.2lb	
<b>Accuracy</b>	50kg±0.2kg; 100kg±0.3kg; 150kg±0.4kg; 200kg±0.5kg; 250kg±0.8kg	
<b>Working Environment</b>	Temperature: 5 C to 40 C    Relative Humidity: ≤90% RH Atmospheric pressure: 86kPa to 106kPa	
<b>Storage Environment</b>	Temperature: -20 C to 60 C    Relative Humidity: ≤90% RH Atmospheric pressure: 50kPa to 106kPa	
<b>Power Source</b>	6V (4 × 1.5V AA alkaline batteries)	
<b>Battery Life</b>	Around 800 measurements (may vary based on usage conditions and network conditions)	
<b>Auto-ON</b>	Sense on technology	
<b>Auto-OFF</b>	About 10 seconds	
<b>Accessories</b>	1. Four AA-size Batteries 2. User Manual	
<b>Mode of Operation</b>	Continuous Operation	
<b>Degree of protection</b>	Type BF applied part (top surface)	
<b>Frequency bands</b>		<b>Transmitted power</b>
GSM900: 880-915MHz/925-960MHz		33dBm+/-2dB
DSC1800: 1710-1785MHz/1805-1880MHz		30dBm+/-2dB
CatM1 & NB IoT 3: 1710-1785 MHz/1805-1880 MHz CatM1 & NB IoT 8: 880-915 MHz/925-960 MHz CatM1 & NB IoT 20: 832-862 MHz/791-821 MHz		23dBm+/-2dB

### ♥ Maintenance

When carrying out usual maintenance, please ensure practice of the following Do's and Don'ts:

- DO use a dry soft cloth to wipe the dust off of this device.
- DO use a wet soft cloth, dipped into water and wrung out, to wipe the dirt off of this device. Then use a dry soft cloth to dry up the device.
- DON'T wash the device with water or immerse it in water.
- DON'T use propellant, abrasive or other chemicals to wipe the dirt off in avoidance of discolor or malfunction.
- DON'T disassemble this device. If you have any problems with this device, please contact customer support.
- Don't perform maintenance while device is in use.

## Warranty

### ♥ Warranty

- Your scale is guaranteed to be free of material and workmanship defects for a period of two (2) years from the date of receipt.
- This warranty does NOT cover damages caused by misuse or abuse, including but not limited to:
  - Failure caused by unauthorized repairs or modifications;
  - Damage caused by shock or drop during transportation;
  - Failure caused by improper operation inconsistent with the instructions stated in this user manual;
  - Malfunction or damage from failure to provide the recommended maintenance;
  - Damage caused by improper use of power supply.

Please contact your program's customer support. Transtek will not be able to directly assist customers.

### **EU regulatory conformance**

Hereby, Zhongshan Transtek Electronics Co., Ltd. declares that this device BS-2001-G1 is in compliance with the essential requirements and other relevant provisions of Directive 2014/53 /EU. The full text of the EU declaration of conformity and most recent information about accessories & software are available at the following internet address: [www.transtekcorp.com](http://www.transtekcorp.com)

### **Manufactured at:**


Zhongshan Transtek Electronics Co.,Ltd  
No.23,Jin'an Road, Minzhong,Zhongshan,528441,Guangdong China

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Medisante Group AG